

## Solutions@Mecmesin Health Care Products

### Specification

Salts Healthcare, a leading manufacturer and supplier of health care products required a low capacity system for testing the integrity of their hydrocolloid products.

### Solution

Salts Healthcare chose the Mecmesin MultiTest 1-*i* because of its accuracy at low loads, ease-of-use and adaptability. They use the MultiTest 1-*i* to conduct quality tests on their hydrocolloid products. Tack testing tests the adhesive qualities of hydrocolloids ensuring that the hydrocolloid products adhere adequately to the skin. The MultiTest 1-*i* is used to test samples every 30 minutes and the results are used by R&D and Production Control to ensure hydrocolloid product integrity.


### System

- MultiTest 1-*i*
- 10N loadcell  $\pm 0.01N$
- Standard accessories

### Testimonial

“The MultiTest 1-*i* has been purchased to test the physical properties of raw materials and finished products. The 10N loadcell is able to measure to two decimal places, which gives a very sensitive and accurate reading, and this test works very well. The MultiTest 1-*i* meets our needs and with time has proved user friendly.”

A Jamesova  
Salts Healthcare



#### Mecmesin Limited

Newton House, Spring Copse Business Park, Slinfold,  
West Sussex, United Kingdom, RH13 0SZ.

sales@mecmesin.com t: +44 (0) 1403 799979 f: +44 (0) 1403 799975 www.mecmesin.com

## medical industry

